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Kirjandus	1. BIOFIRE® FILMARRAY® Torch Kasutusjuhend (2025) 2. BIOFIRE® FILMARRAY® Joint Infection (JI) Panel, ET, REF, RFIT-ASY-0138 3. Proovivõtu käsiraamat 4. Revisiting diagnostics: multiplex infectious disease panels, <i>Clinical Microbiology and Infection</i> 32 (2026) 203-206 5. Tsai SW, Tarabichi S, Vasarhelyi EM, et al. 2025 ICM: Samples for Culture. <i>J Arthroplasty</i>. 2025;41(1 Suppl 1):S277-S284. 6. Esteban J, Patel R, Aguilera-Correa JJ, et al.; Culture Working Group of the Unified PJI Definition Task Force. Optimized use and performance of culture for periprosthetic joint infection diagnosis: a comprehensive literature review. <i>Clin Microbiol Rev</i>. 2025;38(4):e0005425. 7. Martinazzi BJ, Indelli PF, Azboy I, et al. 2025 ICM: Diagnostic Techniques: Molecular Tests. <i>J Arthroplasty</i>. 2025;41(1 Suppl 1):S407-S411. 8. Cashman J, Mortazavi SMJ, Indelli PF, et al. 2025 ICM: Debridement, Antibiotics, and Implant Retention. <i>J Arthroplasty</i>. 2025;41(1 Suppl 1):S297-S325. 9. Indelli PF, Totlis T, Lovreković B, et al. Molecular diagnostics for perioperative microbial identification in periprosthetic joint infection: A scoping review and proposal of a diagnostic flow chart. <i>J Exp Orthop</i>. 2025;12(2):e70263. 10. Ghirardelli S, Scaggiante F, Troi C, et al. Multiplex PCR in septic arthritis and periprosthetic joint infections microorganism identification: Results from the application of a new molecular testing diagnostic algorithm. <i>J Exp Orthop</i>. 2024;11(3):e12097.
Koostaja	Marina Ivanova